


# *Visual Parts and Accessories Listing by Chair Model*

## *Broda Seating*

### **Elite Chairs: Models 785, 35V and 85V Standard Padding**



**N-1**  
Neckrest Cushion  
-Supports Head and Neck.



**SB-1**  
Shoulder Bolster (Pair)  
-Wedges against shoulders to support upper body.



**PMHB**  
Back Pad with terri cloth cover.  
-Additional padding for back.  
-Terri cover is washable.



**SPC**  
Seat Pad with terri cloth cover.  
-Additional padding for seat.  
-Terri cover is washable.



**SP-5U**  
Upper side pad. (Pair)  
-Cushions inside of chair  
around hips.



**SP-5L03**  
Lower side pad (Pair).  
-Covers adjustable arm pins.



**N-6C**  
Calf Pad for Flipdown Footrest.  
-Cushions calves and lower legs.



**N-6S2**  
Sole Pad for Split Flipdown Footrest (Pair).  
-Cushions feet.

**CASTER5''PT**

5'' Plastic Total Lock.

Front caster (used on Elite chairs  
Beginning July 1, 2001)

-Foot brake for safety.

**CASTERT6**

6'' plastic single directional lock


Back caster (used on Elite chairs  
beginning July 1, 2001)

-360° swivel, locks, and tracks.

For previous style Elite casters, see section [Casters-Wheels](#)**Broda Basic Chairs: Models 80V and 30VT Standard Padding****N-1**

Neckrest Cushion.

-Supports Head and Neck.


**SB-1**

Upper Bolster pads (Pair).

-Wedges against shoulders to  
support upper body.**SP-4L**

Lower Side Pad Cover.

-Covers adjustable arm bolts.

**N-6S1**

Sole Pad for Single Flipdown Footrest.

-Cushions feet used on F-FD.

**CASTER5''PT**

5'' plastic total lock

-Foot brake for safety.

(Basic chairs beginning  
July 1, 2001)**CASTER5''PD**

5'' plastic Directional Lock

-Straight rolling caster w/footbrake.


(Basic chairs beginning July1, 2001)

For previous style Basic casters see section [Casters-Wheels](#)

## **Broda Midline Full Recliner**



**N-1**  
Neckrest Cushion.  
-Supports Head and Neck.



**SB-1**  
Upper Bolster pads (Pair).  
-Wedges against shoulders to support upper body




**SPC**  
Seat Pad with terri cloth cover  
-Additional padding for seat  
-Terri cover is washable.



**PMHB**  
-Cushions calves and lower legs  
Terri cover is washable.  
-Additional padding for back



**N-6C**  
Back Pad with terri cloth cover.  
Calf Pad for Flipdown Footrest



**N-6S1**  
Sole pad for Single Flipdown Footrest  
-Cushions feet used on F-FD



**Midline Upper Side Pad**  
(Pair) Cushions inside of chair around hips



**Midline Lower Side Pad**  
(Pair) Covers outer side of chair and armrest adjustment pins



**CASTER5"PT**  
Foot brake for safety.  
5" plastic total lock



**CASTER5"PF**  
5" Rear Fixed Locking Caster.  
Foot brake for safety

**Broda Pedal Chair: Model 48 Standard Padding**



**N-1** Neckrest Cushion.  
-Supports Head and Neck.



**SP-4L** Basic Lower Side Pad (Pair).  
-Covers adjustable arm bolts.



**CASTER5"PF**  
5" Rear Fixed Locking Caster.  
-Foot brake for safety.

**Broda Glider 100-20AL, 100-15AL, 100-10AL:**



**CASTER5"PT**  
5" Front-Locking Caster.  
Freewheeling castor (360°).



**N-1**  
Neckrest cushion  
Supports head and neck

Not Pictured (for Glider)  
**Seat and Back Pad**  
**Side Pads**

For previous style Pedal casters, see section **Casters-Wheels**

***Accessory Cushions: (not applicable to Commode/Shower chairs)***



**N-1 + 2 N-7's**  
Neckrest Cushion (N-1)& two lateral Bolsters (N-7). Supports head and neck.  
Apply to: All models.



**SB-2**  
Extended Shoulder Bolster  
-Part of HSP package; cushions shoulder area.  
Apply to: All except Glider.



**WC-3**  
Short Wing Cover for removable wing (Pair). Part of HSP package. Cushions elbow and hip area of chair  
Apply to: All except Glider.



**WC-4**  
Long wing cover for removable wing (Pair).  
Apply to: All except Glider.



**N-7 Round Bolster**  
-Torso or head support.  
Apply to: All Models.



**N-12 Wedged Bolster**  
-Used in various ways to position resident.  
Apply to: All Models.



**N-6B** Extended Footpad  
-Cushions calves, feet and lower legs.  
Apply to: All except for Pedal and Glider.



**N-9 Legrest Apron**  
-Cushion lower legs.  
Apply to: All models except Pedal and Glider



**N-6BX** Full Extended Footpad  
-Cushions calves, feet and lower legs.  
-part of HSP package  
Apply to: All except for Pedal and Glider.



**HSP** Cushion Package = SB-2, WC-3, ST-4, N-6BX.  
-Extra padding for aggressive behaviors  
Apply to: All except Pedal and Glider.



**SR-1** 1" Side Reduction Pads (Pair).  
Decreases seat width.  
Apply to: All except Glider.



**A-5** Padded Armrest  
Provides arm support and  
Protection for aggressive  
behaviour  
Apply to: All models



**HSP UPGRADED** Cushion  
package = Huntington 's Side Pads,  
Long wing Covers, (not shown in  
picture) Full Extended Footpad, Fixed  
Footrest, Thigh Belt, Fleece Legrest  
Apron  
Apply to: All except Pedal and Glider



**N-1C** Headrest Cushion  
Contoured- **Standard foam**  
insert, contoured for comfort

**N-1CV** Headrest Cushion  
Contoured- **Visco foam** insert,  
contoured for comfort

Apply to: All models



**NS** Neck Support Half  
Circle (Visco Foam)  
Soft Visco foam insert,  
provides comfortable neck  
and head support

Apply to: All models



**SAS** Side Arm Support Covers gap  
between wing and arm  
Apply to: Basic, Pedal and Elite 35V

## Seat Pads (not applicable to Commode/Shower chairs)



**SPC** Seat Pad with removable Cover. Standard with Elite chairs  
Apply to: All



**ESPC** Removable seat cover. Convenient for laundering.  
Apply to: All



**PMHB** Back pad with cover. Back Pad with removable cover.  
Apply to: All



**EHBC** Extra Back Pad Cover. Convenient for laundering.  
Apply to: All



**FWS-SP** Fleece Seat Pad. Soft comfortable surface for skin.  
Apply to: All



**FWS-BP** Fleece Back Pad. Soft comfortable surface for skin.  
Apply to: All



**FWS-FP** Fleece Full Length Pad. Soft comfortable surface for skin.  
Apply to: All



**SPCV** Seat Pad with Removable Vinyl Cover  
**PMHBV** Back Pad with Removable Vinyl Cover  
Apply to: All



**A-FWS** Fleece Armrest Covers (Pair). Soft comfortable surface for skin.  
Apply to: All

## Seat Belts



### **S-2 (Side Buckle top), S-3 (Center Buckle bottom)**

Unpadded Seat Belt 2" wide.

Provides safety; buckle placement allows resident to release buckle.

Apply to: All



### **S-6 Padded Seat Belt 6" wide**

Pads midsection of resident

Apply to: All except Commode/Shower chairs



### **ST-4 Padded Thigh Belt 6" wide.**

Support across the thighs instead of mid-section of resident.

Apply to: All except Commode/Shower chairs

## Trays (not applicable to Commode/Shower chairs)



### **T-8 16", 18", 20" or 22" width**

Slide-on ABS Dinner Tray

Great for activities

Apply to: All



### **T8-C 16", 18" or 20" width**

Clear Dinner Tray

When applied to pedal chair resident can see their feet

Apply to: All



### **T-6L or T-6R Half Tray** Left or

right side available

ABS plastic, helps provide arm support, flips down to one side and is fully removable.

Apply to: Elite, Midline, Glider and 587



**T-10 16", 18" 20" or 22" width**  
Padded T-8 Tray Cover.  
Provides protection and comfort.  
Apply to: All



**T-4 Custom Tray**  
Custom Width Wood Laminate Trays (oversized).  
Made for custom width chairs 14", 24", 26", 28"  
Apply to: All

## **Other Accessories**



**H-4**  
Broda IV Pole: Can be mounted on both sides of chair.  
Apply to: All except Commode/Shower



**PAB BAR**  
PAB Cylinder Length.  
Adjustment Bar.  
Levels seat height of pedal chair.  
Apply to: All except Glider and Commode/Shower



**F-1 (8"), F-2 (12"), F-3 (16")**  
Fixed Footrest.  
Solid footrest with varying lengths.  
Apply to: Basic, Elite, Midline



**F-FD**  
Flipdown Footrest – Single Footplate.  
One solid foot plate that flips up.  
Transfer. Apply to: Basic, Elite, Midline



**F-SD**  
Flipdown Footrest – Split.  
Footplate. Easy front transfer.  
Accommodates feet independently.  
Apply to Basic, Elite, Midline:



**ACS-02**  
Oxygen Tank Holder.  
Holds large and small tanks.  
Apply to: All except Commode/ Shower



**F-SW** Swing Away Footrests (pair)  
Allows resident to rest feet while propelling with mag wheels  
Apply to: Pedal and 587 Wheelchair



**F-SWELV** Elevating Swing Away Footrests (pair) Allows resident to rest and elevate feet while propelling with mag  
Apply to: Pedal



**HEEL LOOPS** (pair) for Swing Away Footrests  
Apply to: Pedal and 587 Wheelchair



**CALF STRAP** for Swing Away Footrests. Helps keep residents' legs in position. **H Strap** also available  
Apply to: Pedal and 587 Wheelchair



**LEGREST PANEL** and **FOOTBOARD**  
Two separate items, shown as one piece (custom) (for Swing Away Footrests) Legrest Panel helps keep residents' legs in position. Footboard attaches to one footplate and flips up for footplate swing away or removal  
Apply to: Pedal and 587 Wheelchair



**LATERAL TORSO SUPPORTS**

Foam with nylon removable cover- Fully adjustable, swing away, comes with mounting hardware.

Apply to: All except Commode/Shower


## Replacement Cushions



**N-1** Neckrest Cushion  
-Supports Head and Neck.  
Apply to: All



**N-1X** Neckrest Cushion Double Thick  
-Supports Head and Neck.  
Apply to: All



**SB-1** Shoulder Bolster (Pair)  
-Wedges against shoulders to support upper body. Apply to: All



**SP-4L**  
Basic Lower Side Pad (Pair).  
Covers arm adjustments bolts.  
Apply to: Basic and Pedal.



**SP-4U**  
Basic Upper Side Pad (Pair).  
Cushioning for hips.  
Apply to: Basic and Pedal.



**SPG-1**  
Glider Side Pads (Pair).  
Cushioning for hips.  
Apply to: Glider (non Auto Locking)




**SP-5L**  
Elite Lower Side Pad (Pair).  
Covers arm adjustments bolts.  
Apply to: All Elite Models.



**SP-5U**  
Elite Upper Side Pad (Pair)  
Cushioning for hips.  
Apply to: All Elite Models.



**N6-C**  
Calf Pad for Flipdown Footrest.  
Cushions calves and lower legs.  
Apply to: All Elite Models.



**N-6S1**  
Sole Pad for Single Flipdown Footrest. Cushions feet.  
Apply to: Basic.



**N-6S2**  
Sole Pad for Split Flipdown Footrest. Cushions feet.  
Apply to: Elite.



**N-66**  
Contracture Sling  
Creates basket area for feet to rest in.  
Apply to: 85V, 35V, 80V and 30VT.

## Replacement Parts



**A-3** Armrest  
For arm support.  
Apply to: All



**W-1R** Removable Wing  
(Specify left or right side).  
Supports upper body  
Apply to: All except Glider and  
Commode/Shower



**W-1** Fixed Wing Frame with mounting  
hardware. Upper body support.  
:



**VSTRAP**  
Vinyl Replacement Straps (with 2 rivets).  
used on all strapped areas of the chair.  
Apply to: All



**A-6** Wooden Armrest  
Apply to: All except Commode/Shoer

## Gas Springs



**PAB-4**  
Gas Spring 300N  
(short model).  
Apply to: 785, Pedal,  
Commode/Shower




**PAB-5**  
Gas Spring 300N  
(standard model).  
Apply to: All except  
Pedal and Glider



**Cable Release**  
Remote cable release  
assembly for PAB (handle,  
cable, actuator).  
Apply to: All except Glider.



**PRT-CABHAND**  
Handle and cable for PAB release Assembly  
Apply to: All except Glider.



**Cable Actuator**  
Actuator for PAB release assembly.  
Apply to: All except for Glider

## Casters - Wheels



**CASTER5"TL**  
5" Total Lock Caster.  
Apply to: Basic prior to  
July 1, 2001



**CASTER5"FX**  
5" Fixed Locking Caster.  
Apply to: rear of Pedal  
prior to Nov 21, 2001.



**CASTER5"NL**  
Non-Locking Caster.  
Apply to: front of Pedal  
prior to Nov 21, 2001



**CASTER5"DL**  
5" Directional Locking Caster.  
Apply to: Basic  
Prior to July 1, 2001  
(Note: Has a blue foot pedal)



**CASTER4"SL**  
4" White Total Lock Front Caster.  
Apply to: Elite Models  
May 28, 1998 to May 15, 2000



**CASTER5"SL**  
5" (Elite line) Total Lock Front Caster.  
Apply to: Elite Models.  
beginning May 15, 2001  
to July 1, 2001



**CASTER6"SL**  
6" Single Directional Locking Combi Caster  
Apply to: Elite Models  
Since May 28, 1998  
to July 1, 2001



**WHEEL24**  
24" Replacement Wheelchair Plastic handrim.  
Apply to: All except Glider Chair.



**WHEEL20**  
20" Replacement Wheelchair Aluminum handrim.  
Apply to: All except Glider Chair



**CASTER5"PF**  
5" rear plastic fixed locking caster  
Apply to: Pedal beginning Nov 21, 2001



**MDLC8"**  
8" Multi-Directional Rear Caster.  
Apply to: Elite Models  
prior to 1997.



**CASTER5"PT**  
5" plastic total lock. Commode shower chair,  
(Basic and Elite chairs beginning July 1, 2001, Pedal beginning Nov 21, 2001)  
Front Caster



**CASTER5PD**  
5" plastic directional lock Commode Shower chair,  
(Basic chairs beginning July 1, 2001)  
Rear Caster



**CASTERT6**  
6" plastic single directional lock. Apply to: Elite chairs  
beginning July 1, 2001  
Rear Caster

**Tray Straps**



**T-3 1"**  
Replacement Strap for T-8.



**T-3 2"**  
Replacement Back Strap for T-2, T-4.

**Replacement Cushions (Before 1997)**



**SP-3**  
Elite Side Pads for 94 to 97 Models (Pair).  
Apply to: Elite Models.



**N-6**  
Footpad for F-1, 2, 3 Footrest.  
Apply to: Elite and Basic Models.



**N-9**  
Legrest Apron.  
Apply to: Elite and Basic Models.

**Anti Tip Stabilizers**



**FCP-RSTABILEL**  
Rear Stabilizers  
Apply to: Elite Models



**FCP-STABILIZER**  
Lateral Stabilizer  
Apply to: Elite, Basic and Pedal



**FCP-REARSTABIL**  
rear stabilizers  
**FCP-FTMSTABIL**  
front stabilizers (look similar to rear stabilizers)  
Apply front to: Elite, Basic and Pedal  
Annlv rear to: Basic and Pedal

**Commode/Shower Chair: Model CS 385**



**CS-PD-SEAT Padded Commode Seat**  
Urethane foam with durable 4 way stretch vinyl cover, securely attaches to chair




**Extended Arms**  
Interchangeable with standard arms, increase seat width up to 25"



**CS PAN Commode Pan**  
Plastic commode pan for **padded** commode seat. Slides on with Wire pan holder



**CS BP BSKT2** Wire commode pan holder. For use only with **padded** commode seat for CS PAN



**CS HDRST**  
Adjustable headrest to support neck and head




**CASTER5"PT**  
5" plastic total lock



**CASTER5PD**  
5" plastic directional lock for tracking



**S-2 (Side buckle, top), S-3 (Center buckle, bottom)**  
Unpadded Seat belt 2" wide. Provides safety; side buckle placement allows resident to release buckle.



**CS-SL-PD**

Side Lateral Supports for CS 385.  
Can also be used on seat



**HDRST-PD**

Adjustable padded headrest to support  
neck and head

# BRODA Seating

April 2004

(Division of BRODA ENTERPRISES INC.)  
 385 Phillip Street., Waterloo, ON N2L 5R8  
 PO Box 33046, G.P.O. Station, Detroit, MI 48232  
 Phone: (519) 746-8080 Toll Free 1-800-668-0637  
 Fax: (519) 746-8616

## WARRANTY GENERAL INFORMATION

BRODA provides a Three Year Warranty on the chair frame, and a One Year Warranty on all other components subject to the following conditions:

No warranty is provided on seat pads or cloth covers.

The chair frames are guaranteed for three years against structural defects or failure. All other parts (except seat pads and cloth covers) including but not limited to strapping, cushions, gas springs and attachments, casters, wheels, brakes and armrests are guaranteed for one year against defects in materials and workmanship based on normal institutional use. The guarantee does not cover malicious or deliberate damage or damage from misuse.

Modifications to Broda products or the use of non-Broda supplied parts voids the warranty. This warranty does not cover shipping damage (see below).

BRODA will provide new or refurbished parts for installation by the owner at no cost following confirmation by the local BRODA Representative or the BRODA Head Office Customer Service Representative. On request, defective parts must be returned to the factory within thirty days of receipt of the replacement parts by the owner. If the defective parts are not returned to BRODA on request, the owner will bear the cost of the replacement parts on invoice from BRODA.

Warranty does not include on-site labor for the installation of warranty parts or warranty repairs. The owner may return to BRODA products for warranty replacement or repair by shipping items prepaid and insured to the factory. Warranty completed at the factory includes both materials and labor. Parts to be repaired or replaced at the discretion of BRODA. All returns to the factory require prior authorization from BRODA.

BRODA retains the right to make design and application changes without notice. All orders will be filled with BRODA's current models unless otherwise specified by the purchaser.

BRODA chairs are designed for patient mobility, positioning, and comfort in specialty seating, however, the application of BRODA products shall remain the responsibility of the purchaser or user.

This warranty is not transferable.

## RETURNS

BRODA Seating will not accept any returns without a prior **Returned Goods Authorization Number**. Please contact our Head Office Customer Service Representative at 1-800-668-0637 for assistance. Returns must be insured when shipped.

## DAMAGED FREIGHT

### NOTIFY THE CARRIER OF ANY DAMAGE IMMEDIATELY

It is the responsibility of the person receiving the goods to examine cartons and goods before accepting receipt. Note all damages on the bill of lading and file a claim if necessary. Notify the carrier of any concealed damaged within a maximum of 48 hours.

BRODA insures all products for intransit damage, failure to notify the carrier of intransit damage voids both the insurance and the BRODA warranty. If you require assistance, please contact our Head Office at 1-800-668-0637.

### SALES TAX

Most Broda products are G.S.T. Zero Rated and Exempt from Canadian Provincial Sales Tax. Purchasers maybe required to check with their Provincial or State Tax Office for purchaser tax payment.