

# BRODA<sup>®</sup> SEATING

*Where Comfort is Essential*



Ideal for nursing homes and various hospital settings including Acute Care, Transitional Care, Day Surgery, Emergency Wards, Orthopedics

## STANDARD FEATURES

- 20" or 22" seat width
- 18.5" seat depth
- 20.5" seat height
- 32" back length
- 6" multi-directional locking rear casters
- 5" total lock front casters
- 16 gauge tubular steel frame
- Infinitely adjustable 27° seat tilt, 90° recline
- Removable fully padded seat and back
- Broda Comfort Tension Seating™
- Infection control friendly full padding package
- Upper lateral supports (wings)
- Removable, swing away sides
- Height adjustable arms
- Elevating and length adjustable legrest

Options include Dinner Tray, Half Tray, Padded Tray Cover, IV Pole, O2 Holder, Mag Wheels, Padded Armrests

## Acute Care Chair

### Ultimate Comfort and Support

Designed to Offer Superior Comfort, Functionality and "Infection Control- Friendly" Cleaning- All in One Chair

- Broda's Comfort Tension Seating™ aids in pressure redistribution, creating *optimal patient comfort*
- Infinitely adjustable tilt, recline and legrest allow for *easy positioning for multiple patients*
- Minimized stitching and removable padding allows for easy, thorough cleaning, thus *accommodating infection control concerns*
- Durable padding, parts and 16 gauge tubular steel frame *reduce chair maintenance and replacement cycles*



### Up to 27° of tilt in space:

- Helps prevent falls
- Enhances patient comfort
- Helps position patients who have minimal upper body control

### Up to 90° of back recline (Trendelenburg capable)

- Assists with patient transfers
- Patients can sleep in the chair

Broda Seating, 560 Bingemans Centre Drive, Kitchener, Ontario, Canada N2B 3X9

Tel: 519-746-8080 Fax: 519-746-8616 Toll Free US/Can. 1-800-668-0637

www.brodaseating.com • sales@brodaseating.com

*A Multi Use Chair*



**Broda Comfort Tension Seating**

*The key to comfort is in the strapping:*

- Reduces pressure, thus helping to reduce skin breakdown
- Maximizes Comfort
- Will not bottom out



Strapping returns to its original shape within 10 seconds after weight is removed

*Patients can receive treatments in the chair as well as using it for sleeping, eating and socializing.*

*Caregiver and Infection Control Friendly!*

*The Acute Care Chair is covered to allow for quick and thorough cleaning.*

- Stitching is minimized to help prevent fluids and bacteria from entering
- Vinyl has an Antimicrobial, Antibacterial and Anti-stain finish.
- Seat, back, sides, wings and padding are removable for cleaning

*For everyday cleaning, use a straight application of a concentrated cleaner or a hospital grade cleaner. Follow diluting instructions on hospital grade cleaner label. Wipe with a clean cloth. Remove parts and padding as necessary for thorough cleaning*



Swing away removable sides allow for easy, thorough cleaning as well as ease of side transfers



Legrest folds under seat to assist patients getting in and out of chair



Legrest is length adjustable to accommodate various patients

**Place Business Card Here**